technique

0123

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Intended use

SinCrest is a kit for the maxillary sinus lift with crestal approach.

The SinCrest system consists of a kit including electric hand-piece drills, depth stops at different levels and a manual osteotome.

The drills and stops allow obtaining a guiding hole in the alveolar bone in the closest possible position to the membrane of the maxillary sinus floor.

The manual osteotome was designed to obtain the controlled fracture of the bone floor near the maxillary sinus membrane without damaging the membrane itself. Thus dental implants can be inserted exceeding the initial bone length.

- **SinCrest Kit is not sterile when supplied. It must be washed and sterilised every time before using it.**
- **The device must be used exclusively by skilled medical staff only.**
- **The surgeon must assess whether the patient is suitable for osteotomy surgery and choose the correct surgical sequence.**
- **When handling the device, always use sterile gloves.**

- **Never use SinCrest without making sure that all its components have been assembled correctly (see fig. 8)**
- **Before using SinCrest, press the probe onto a sterile field and check the correct recovery by means of the spring. (see fig. 9)**
- **Make sure that a dental floss is inserted onto the RX Pins to avoid the risk of accidental ingestion by the patient.**
- **An incorrect speed (RPM) of the hand-piece may overheat the bone tissue leading to bone necrosis.**
- **META does not shoulder any responsibility for the incorrect use of the product.**

Contraindications

- Do not use the SinCrest technique in case of residual bone shorter than 5mm and taller than 11mm.
- The Ref 4344 Probe Drill Ø 1.2 is not designed to cut the bone laterally and must be used only to drill in depth. An incorrect use may damage the drill.
- We recommend not to use the SinCrest in case of reduced stability during the screwing stage in the bone hole drilled.

Usage procedure

Precautions

Correct X-rays, direct palpation and inspections of the operating site are necessary to identify the anatomy of the available bone. The most suitable surgical sequence is to be assessed according to the available bone and the size of the dental implant.

During the use of the drills to prepare the bone hole, ensure an abundant external lavage with sterile saline solution; a good lavage may also be ensured by means of a sterile syringe.

Surgical procedure

After having proceeded to the full depth incision and elevation of the soft tissues, the site may be prepared to use the selected SinCrest device.

To handle the drills, we suggest to use the accessory called "Drill Gripper", which allows preventing the direct contact with the operator's hands. From the side marked as N. 4, the Ref 4381-SinCrest drill Ø 4.0 and Ref 4340-SinCrest drill HB Ø 4.2 may be withdrawn, whereas the other side contains all the other drills of the SinCrest Kit (see figure n. 1)



Fig. 1

The laser marks reported on the Probe drill Ø 1.2 and the SinCrest drills Ø 3.0, 3.2, 4.0 and 4.2 indicate the operating depth at 4-5-6-7-8-9-10 and 11mm. (see figure n. 2)

Surgical sequence

Ref 4343-Locator drill: use the drill to identify the correct position of the implant and remove the cortical bone tissue at a 3.5mm depth as shown by the laser mark.

Ref 4344-Probe drill Ø 1.2: to be used after having positioned the depth stop selected according to the bone availability (see paragraph on Stop positioning). Since the osteotomy achieved defines the final position and angle of the implant, it must be carried out with a view to the full prosthetic plan. Use the drill with back and forth movements until the desired depth is achieved. Only the head of the device is sharp. No lateral cutting strength can be implemented to modify the osteotomy angle: it must be used only to drill in depth.

Ref 4354-RX pin Ø 1.2: use the pins to check the correct angle and depth achieved through an x-ray and/or visual examination.

Ref 4379-Sin Probe Ø 1.2: in addition to the tactual check, this device allows a metric assessment of the operating depth.

Ref 4380-Guide drill Ø 3.0: use the drill to create a coronal recess in the cortical layer to achieve a 2 mm depth as marked on the drill.


Ref 4337-SinCrest drill Ø 3.0: the use of this drill must always be preceded by the Guide drill. It may be used after having positioned the depth stop selected

according to the available bone (see paragraph on Stop positioning). Prepare the site until the pre-determined depth is achieved.


Ref 4356-RX pin Ø 3.0: confirms the correct angle and depth achieved during the previous osteotomy stages.

Ref 4333-SinCrest Ø 3.0:


Step 0 Gradually screw the device inside the dedicated guide hole. The device has achieved the depth obtained by means of the drills when the black mark emerges minimally during the progress of the probe as compared to the grip (see figure n. 3).


Step 1 Turn the handle anti-clockwise by 1/2 revolution and at the same time implement an axial pressure, start the osteotomy with "back and forth" turning movements by 1/2 turn until the distance (0.5mm) between the tooth and the stroke is zero (see figure n. 4 - 5 A).

By manually implementing pressure on the probe (see figure n. 5 B), constantly check the residual strength of the maxillary sinus floor.

Step 3 Screw by 1/2 turn (each turn equals a 0.5 mm progress) after the tooth is engaged in the stroke (see figure n. 6) and repeat the same steps until the fracture of the maxillary sinus floor detected by the probe advance in the apex direction is reached.

Ref 4338-SinCrest drill HB Ø 3.2: this drill must be used in case of hard bone encountered during the screwing stage of the SinCrest Ø 3.0 device. This drill does not require positioning a stop since it must simply widen the hole without working in depth. Prepare the site until the desired depth has been achieved. The osteotomy may be resumed by means of the SinCrest Ø 3.0 device. Proceed to insert the bio-material or the dental implant directly. Restore the position of the flap and suture.

If bone is available, the surgical procedure may be continued as illustrated in the following steps so as to position a dental implant with a wide diameter.

Ref 4381-SinCrest drill Ø 4.0: the use of this drill must always be preceded by the SinCrest drill Ø 3.0. This drill does not require a stop since it must simply widen the hole without working in depth. Prepare the site until the desired depth is achieved by means of the SinCrest drill Ø 3.0.

Ref 4357-RX pin Ø 4.0: confirms the correct angle and depth achieved during osteotomy.

Ref 4334-SinCrest Ø 4.0:

Step 0 Gradually screw the device inside the dedicated guide hole. The device has achieved the depth obtained by means of the drills when the black mark emerges during the progress of the probe as compared to the grip (see figure n. 3).

Step 1 Turn the handle anti-clockwise by 1/2 revolution and at the same time implement an axial pressure, start the osteotomy with "back and forth" turning movements by 1/2 turn until the distance (0.5mm) between the tooth and the stroke is zero (see figure n. 4 - 5 A).

By manually implementing pressure on the probe (see figure n. 5 B), constantly check the residual strength of the maxillary sinus floor.



Step 3 Screw by 1/2 turn (each turn equals a 0.5 mm progress) after the tooth is engaged in the stroke (see figure n. 6) and repeat the same steps until the fracture of the maxillary sinus floor detected by the probe advance in the apex direction is reached.

Ref 4340-SinCrest drill HB Ø 4.2: this drill must be used in case of hard bone encountered during the screwing stage of the SinCrest Ø 4.0 device. This drill does not require positioning a stop since it must simply widen the hole without working in depth. Prepare the site until the desired depth is achieved. The osteotomy may be resumed by means of the SinCrest Ø 4.0 device. Proceed to insert the biomaterial or the dental implant directly. Restore the position of the flap and suture.

Stop positioning

8 stops are available for different depths (4-5-6-7-8-9-10 and 11mm) with laser marks for easier identification. The stops are organised in sequence in the dedicated housings in the tray. Their positioning on the Probe Drill Ø 1.2 and SinCrest drill Ø 3.0 takes place directly by simply inserting the drill (see figure n. 7) mounted on the contra-angle hand-piece, directly in the desired depth stop. The stop already used may be temporarily positioned in the stop housing, thus avoiding the direct contact with the operator's hands.

To avoid the accidental fall of the stops, slide the cover into the lock position.

Instructions to use the sterilisable surgical devices

Washing and disinfecting


SinCrest Kit is not sterile when supplied. It must be washed and sterilised each time before it is used.

Disinfect the devices by means of the dedicated germicide detergents immediately after the surgical procedure to prevent cross contamination.

To wash the SinCrest Ø 3.0 and Ø 4.0 device accurately, use the Screwdriver accessory contained in the tray which allows unscrewing the top container and removing the 2 internal components (see figure n. 8). The inner cavity of the device must be washed by means of a brush.

Wash the devices accurately with flowing water to remove the residual tissue with the help of a soft brush and check that they are not damaged.

For ultrasound cleaning, position the surgical devices inside a container separately from other materials so that no damage nor contamination takes place.

Sterilisation and storage

Before proceeding with the sterilisation stage, correctly re-assemble the components and make sure that the spring is repositioned as shown in the figure (greater diameter facing the upper container).

Use the dedicated SinCrest Tray to arrange the surgical devices. Disinfect and wash the Tray each time immediately after it has been used. Make sure that no residues remain in the drill and tool housings; in certain cases, disassembling the tools is recommended for better cleaning. When the Tray is washed in dedicated tool-washers, place the tray diagonally for a better elimination of the condensate steam. The use of incorrect chemical detergents may damage the tray. Seal the tray individually in a suitable sterile bag.

Sterilisation should take place inside a steam steriliser at 134°C / 18 min. / 2,150 bars, always using validated sterilisation cycles according to the laws in force. The maintenance of sterilisation is ensured by the whole bag: always check that it is not damaged before using it.

You can avoid mistakes and ensure the correct operation of our product only if you read the present document.

The present operating instructions and additional information on the SinCrest technique are available on our web site: www.metahosp.com

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List of sterilisation and disinfection devices

SinCrest Kit Ø 3.0 Ref 4417	SinCrest Kit Ø 3.0/4.0 Ref 4418	Product	Ref	Description	Reusable	Steam sterilisation	Disassembling recommended	Material
●	●		4343	LOCATOR DRILL	Yes ^{1,2}	Yes		Stainless Steel
●	●		4344	PROBE DRILL Ø 1.2	Yes ^{1,2,3}	Yes		Stainless Steel
●	●		4380	GUIDE DRILL Ø 3.0	Yes ^{1,2}	Yes		Stainless Steel
●	●		4337	SINCREST DRILL Ø 3.0	Yes ^{1,2,3}	Yes		Stainless Steel
	●		4381	SINCREST DRILL Ø 4.0	Yes ^{1,2}	Yes		Stainless Steel
●	●		4338	SINCREST DRILL HB Ø 3.2	Yes ^{1,2}	Yes		Stainless Steel
	●		4340	SINCREST DRILL HB Ø 4.2	Yes ^{1,2}	Yes		Stainless Steel
●	●		4354	RX PIN Ø 1.2	Yes	Yes		Ti6Al4V
●	●		4356	RX PIN Ø 3	Yes	Yes		Ti6Al4V
	●		4357	RX PIN Ø 4	Yes	Yes		Ti6Al4V
●	●		4379	SIN PROBE Ø1.2	Yes	Yes		Ti6Al4V
●	●		4333	SINCREST Ø 3.0	Yes ^{1,4}	Yes	Yes	Ti6Al4V/ Stainless Steel/ Plastic
	●		4334	SINCREST Ø 4.0	Yes ^{1,4}	Yes	Yes	Ti6Al4V/ Stainless Steel/ Plastic
●	●		4339	DRILL GRIPPER	Yes	Yes		Plastic
●	●		4341	SCREWDRIVER	Yes	Yes		Plastic
●	●		4406 4407 4408 4409 4412 4413	DRILL STOP 4-5-6-7-8-9-10-11mm	Yes	Yes		Titanium
●	●		4342	SINCREST TRAY	Yes	Yes	Yes	Plastic/ Stainless Steel

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