Zmatrix™ - collagen membrane

DESCRIPTION

Zmatrix™ is a biodegradable barrier membrane for use in Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR) procedures. In these procedures, Zmatrix™ can be used as a resorbable membrane for the containment of bone grafts in bone repair procedures during treatment of periodontal bone defects (GTR) and for bone augmentation procedures (GBR). Zmatrix™ is a highly purified membrane produced from porcine collagen using controlled and standardized manufacturing processes. Starting materials for the production of Zmatrix™ are harvested in EU certified facilities under strict veterinary controls from veterinary certified animals that are declared fit for human consumption. Zmatrix™ is packaged in double blisters and sterilized by ethylene oxide gas treatment.

PROPERTIES

Zmatrix™ is composed of a network of highly purified porcine collagen fibers intermingled with highly purified porcine elastin fibers. The membrane is designed to protect the bone defect space in order to establish an environment favorable for bone regeneration by preventing the migration of undesired cells from the surrounding soft tissue and allowing the ingrowth of osteogenic cells. The fiber network of Zmatrix™ provides a physical barrier to contain the bone graft material at the defect site and to exclude in-growth of surrounding tissue for a period of time long enough to allow bone regeneration to take place. Zmatrix™ is biodegradable and a second surgical procedure to remove the membrane is not needed. The biocompatibility of Zmatrix™ has been demonstrated through testing and is due to the utilization of effective purification and processing methods during manufacturing. Experimental studies have shown that Zmatrix™ can be considered safe regarding potential porcine viral transmission. Minor variations in thickness and a slightly irregular collagen structure can be occasionally seen with dry Zmatrix™. This is due to the natural origin of the membrane and does not influence its clinical performance. Zmatrix™ is hydrophilic and retains its structural integrity also in the wet state. Application of Zmatrix™ should involve the use of space-maintaining bone graft material, whenever possible. There is no preferred orientation of Zmatrix™ during application.

INDICATIONS FOR USE

Zmatrix™ is intended for use during the process of Guided Bone Regeneration (GBR) and Guided Tissue Regeneration (GTR) as a biodegradable barrier membrane for:

- Ridge augmentation for later implant insertion
- Augmentation around implants placed in immediate/delayed extraction sockets
- Localized ridge augmentation
- Alveolar ridge preservation/reconstruction
- Osseous fill around implants in peri-implantitis bone defects
- Over the window in lateral window sinus elevation procedures
- Intra-bony defects around teeth
- Treatment of recession defects, together with coronally positioned flap
- In furcation defects in multi-rooted teeth.

INSTRUCTIONS FOR USE

The principles of sterile handling and applicable patient medication must be followed.

Directions for use in GBR/GTR procedures:

- The bone defect is exposed surgically (mucoperiosteal flap) and the soft tissues are removed using standard surgical procedures.
- In case of GTR procedures the surgical treatment should be preceded by a conditioning with cause-related anti-infective therapy to eradicate any remaining bacterial infection and a thorough training of the patient in oral hygiene. During treatment the root surface should be debrided and planed carefully.
- Before bone augmentation around implants in peri-implantitis bone defects can be performed, appropriate debridement and implant surface disinfection should be achieved.
- The bone defect is filled with bone graft material (e.g. autogenous bone, xenogeneic bone, synthetic bone substitutes, or mixtures) which may be mixed with blood, bone marrow, cell fractions or other substances to enhance bone formation.
- A sterile template (e.g. a sterile aluminum foil) may be used to determine the required size and shape of Zmatrix™.
- Zmatrix™ is trimmed to the required size with sterile scissors. The final trimmed Zmatrix™ should overlap the defect walls by at least 2 mm to avoid a lateral in-growth of gingival tissue and to assure a complete bone contact of the membrane.
- The trimmed Zmatrix™ is placed over the defect and its contact to the bone surface and bone graft material can be supported by applying mild pressure until Zmatrix™ adheres completely. A complete wetting of Zmatrix™ with blood or exudate is recommended. Zmatrix™ has excellent handling properties when dry or wet and it can be repositioned if necessary.
- The mechanical properties of Zmatrix™ allow a fixation of the membrane if required to avoid a displacement.
- Finally, the initially generated soft tissue flap (mucoperiosteal flap) is closed completely over the membrane by tensionless sutures in order not to compromise blood supply in the defect area. If the wound cannot be closed completely an accelerated degradation of the membrane and a bacterial infection of the defect area may occur.
- In case of GTR procedures it is important to adapt Zmatrix™ closely to the respective tooth optionally with additional sutures to avoid the formation of excessive junctional epithelium.

Post-operative care:

- The patient should be monitored closely.
- If Zmatrix™ becomes exposed, the dehiscence usually heals by itself within several weeks. Membrane removal is usually not necessary. However, to minimize bacterial contamination rinsing with bactericidal solutions is recommended.
- In the event that early removal of Zmatrix™ is necessary, the tissues adjacent to the membrane should be anesthetized with a local anesthetic. An incision should then be made immediately adjacent to the residual membrane. Following careful reflection of the surrounding tissue, the remaining portion of the membrane can be excised and the area curetted to remove any inflamed or infected tissue.

- To allow for undisturbed bone regeneration underneath Zmatrix™ surgical reentry should not be done before 4 to 6 months postoperatively.

LIMITATIONS OF USE

Contraindications:

- The use of Zmatrix™ without space maintaining bone grafting material is not indicated.
- Zmatrix™ should not be used in infected areas. An active infection of the implant site may lead to an increased biodegradation of Zmatrix™ which may lead to an early loss of its barrier function.
- Adverse reactions during the use of porcine-derived collagen membranes in GBR/GTR procedures have not been observed. However, allergic reactions cannot be totally excluded. The use of Zmatrix™ in patients with known sensitivity to porcine-derived materials and collagen is not recommended.

Undesirable side effects:

- In case of soft tissue wound dehiscence, an increase in biodegradation of Zmatrix™ may occur, which could result in an earlier loss of the barrier function.
- Other possible side effects are not specifically related to Zmatrix™ itself but may occur with any surgery, like infection, wound swelling, bleeding, local inflammation, bone loss, or pain.

Precautions:

- In patients with constrained healing capacity (e.g. due to metabolic disorders, cancer treatments, heavy smoking) Zmatrix™ should be used with special caution and closer patient monitoring during post-operative care.
- Clinical data on the use of Zmatrix™ during pregnancy or lactation is currently not available. As a safety precaution, it is not recommended to use Zmatrix™ in women who are pregnant or breast feeding.

WARNINGS AND OTHER INFORMATION

Zmatrix™ should not be used when the sterile packaging is damaged.

Zmatrix™ should not be resterilized and is for single use only. A resterilization will lead to a deterioration of product properties and may result in unwanted biological reactions due to an unsterile product or to a treatment failure due to an early loss of the membrane's barrier function.

STORAGE AND HANDLING

Zmatrix™ should be stored at room temperature ($\leq 25^{\circ}\text{C}$ / $\leq 77^{\circ}\text{F}$) and in a dry place. Zmatrix™ should be handled using sterile gloves and/or with sterile instruments.

Do not use after the use by date indicated on the packaging.

The manufacturer is not liable for an unintended or incorrect use of Zmatrix™.

PRESENTATION FORM

One unit of Zmatrix™ is provided per double blister package.


Caution

Rx Only

U.S. Federal law restricts this device to sale by or on the order of a dentist or physician.

Date of information: 05/2023 Rev 0523

Glossary of Labeling Symbols

SYMBOL	DESCRIPTION
	Manufacturer
	Use by date
	Lot No.
	Article No.
	Sterilized using ethylene oxide
	Do not resterilize
	Do not use when sterile packaging is damaged
	Upper temperature limit
	Single use
	See Instruction For Use
	Caution. Consult post-operative care and limitations of use.

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In Canada, Manufactured for:

Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424 USA
osteogenics.com
+1-806-796-1923

Distributed by:

Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424 USA
osteogenics.com
+1-806-796-1923

Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79, Lubbock, TX 79424 USA
+1-806-796-1923 | osteogenics.com

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Matricel GmbH