Allotech, LLC - enCore Allograft Package Insert

DESCRIPTION / USE: In accordance with FDA Article 21 CFR Part 1271, this package contains donated human cells, tissues, and cellular and tissue-based products (HCT/Ps). Human allograft bone may be used in a variety of dental reconstructive, regenerative, or periodontal procedures. Tissue is processed and preserved by a variety of techniques, supplied in a range of sizes for surgical use by licensed clinicians and is intended for homologous use only. The tissue is donated human tissue. All tissue is processed and packaged using aseptic technique, freeze dried, and terminally sterilized.

CONTRAINDICATIONS: The presence of infection at the transplantation site is a contraindication for use of musculoskeletal allograft.

WARNINGS AND PRECAUTIONS: Human tissue has the potential to transmit infectious agents. Donor screening, processing treatments and laboratory testing follow stringent specifications to reduce the risk of infectious agent transmission. Do not use if the expiration date has been exceeded or if there is evidence of defects in package or label integrity. Do not sterilize or re-sterilize. It is the responsibility of the clinician to maintain tissue for transplantation according to recommended storage conditions. Do not use if tissue has not been stored according to the recommended STORAGE instructions. This tissue is intended for use by licensed clinicians.

Allografts are processed using some or all of the following agents: povidone iodine, physiological buffers, alcohols and saline solutions. Traces of these reagents may remain on the tissue allografts.

Although this tissue has been tested and screened for relevant communicable diseases and disease agents, and processed under aseptic conditions, application and use of any allograft tissue may potentially have negative outcomes. Occurrence of complications at the affected site may transpire post-treatment, without early warning signs. These include, but are not limited to: 1) transmission of communicable diseases; 2) transmission of infectious disease agents; and 3) immune rejection of, and/or allergic reaction to the HCT/P. Any adverse outcomes potentially attributable to the HCT/P must be reported promptly to Allotech, LLC.

Caution should be taken when administering this product to immunocompromised individuals, such as patients suffering from HIV or other highly immunocompromised conditions. Although Allotech, LLC has taken great measures to ensure the safety of our allograft products, current technologies cannot preclude the transmission of certain diseases known or unknown. Therefore, Allotech, LLC can make no claims concerning the biological properties and safety of allograft tissue.

The end-user shall be held responsible for determining the appropriate application and usage of this product.

DONOR ELIGIBILITY & SEROLOGICAL TESTING: Donor eligibility (screening and testing) is performed in accordance with AATB standards and FDA regulations. Donor screening includes assessment of the medical and social history as well as physician assessment of the donor to assure that no conditions exist that may make the tissue unacceptable for transplantation. Donor determined to be eligible, based on the results of screening and testing. Screening includes a review of medical and social history, available hospital

records, infectious disease testing, autopsy reports (if performed), and physical examination of the Donor, at the time of recovery. The Donor has been tested using FDA licensed, approved, or cleared donor screening test kits and was found negative or non-reactive for:

- Human Immunodeficiency Virus Types 1 and 2 Antibody (anti-HIV-1/anti-HIV-2)
- Hepatitis B Surface Antigen (HBsAg)
- Hepatitis B Core Antibody – Total (anti-HBc)
- Hepatitis C Virus Antibody (anti-HCV)
- Human Immunodeficiency Virus 1, Hepatitis B Virus and Hepatitis C Virus Nucleic Acid Test(s) (HIV 1/HBV NAT/HCV NAT)
- Syphilis Rapid Plasma Reagin or Treponemal Specific Assay

Additional tests, including but not limited to HTLV I/II, CMV or West Nile Virus, may have been performed and, if performed, were found to be acceptable for transplantation. Communicable disease testing has been performed by a laboratory registered with the FDA to perform donor testing in accordance with the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR part 493, or that has met equivalent requirements as determined by the Centers for Medicare and Medicaid Services (CMS).

MICROBIAL TESTING AND TERMINAL STERILIZATION: Prior to packaging, tissue is subjected to microbiological testing at recovery and in the course of processing and must be free of specific aerobic / anaerobic microorganisms and fungal contaminants whose presence would preclude tissue from processing or transplantation. Following packaging, tissue is terminally sterilized by e-beam irradiation to a minimum sterility assurance level (SAL) of 10⁻⁶.

POTENTIAL COMPLICATIONS / ADVERSE REACTIONS: Inherent uncertainty exists in medical and social histories and laboratory testing which may not detect known or unknown pathogens. Therefore, the following complications may occur with tissue transplantation: Loss of function or integrity of transplanted tissue with resorption, fragmentation, disintegration, and associated loss of continuity, displacement, bending or fracture. Immune response to transplanted tissue. Transmission of known pathogens including Hepatitis B or C, Human T-Cell Leukemia / Lymphotropic Virus, Human Immunodeficiency Virus 1 & 2, Syphilis, or bacteria. Transmission or causation of diseases of unknown etiology and characteristics.

BEFORE USE: Examine Allograft Packaging. Do Not Use This Allograft If: Any of the packaging elements appear to be missing, tampered with or damaged. The product label is severely damaged, illegible, or missing. The expirations date shown on the package label has passed.

GRAFT STORAGE: It is the responsibility of the Tissue Dispensing Service, Tissue Distribution Intermediary, and/or End-User clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to further distribution or transplant. Store freeze dried grafts at 15-30°C (59-86°F) in a clean, dry location. Allograft must be stored with proper precautions to prevent contamination.

DIRECTIONS FOR HANDLING AND PREPARATION: All preparation should be performed using aseptic technique. Carefully open the outer tray of the single blister and aseptically remove the sterile inner container holding the ground bone matrix. The sterile ground bone matrix can then be removed from the sterile inner container for usage during the surgical procedure. Once the

packaging has been opened, the tissue must either be transplanted or discarded. The tissue is intended for use in one patient, on a single occasion only.

Graft Type	Graft Storage	Graft Preparation
Ground Bone	Store freeze dried grafts at 15-30°C (59-86°F) in a clean, dry location.	To reconstitute, place graft in sterile basin and cover with sterile isotonic solution for 15-30 minutes, depending on the size of the graft

TISSUE TRACKING: The FDA requires patient records to be properly maintained by storing the allograft LOT NUMBER for purposes of tracking the allograft from the donor to the recipient. Please go to our website, <https://osteogenics.com/fda-tissue-tracking>, and submit tissue tracking information.

LABELING SYMBOLS: Symbols may be used on package labeling for easy identification.

 Expiration Date	 Single use only	 Lot Number
 Attention, see instruction for use	 Method of Sterilization Using Irradiation	
 15°C ↓ 30°C	Temperature limit 15-30°C (59-86°F)	 Catalog Number

Red indicator label indicates exposure to irradiation. The tissue's expiration date is expressed in the format YYYY-MM-DD.

Allotech, LLC processed, stored, and distributed all tissue in accordance with the standards of the American Association of Tissue Banks (AATB), the FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), and applicable State regulations. Osteogenics Biomedical stored and distributed all tissue in accordance with the FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), and applicable State regulations.

Responsible Tissue Banks:
Source Establishment, Determines Donor Eligibility & Processing: Surgenex, LLC, 15444 N. 76th Street Suite #110 Scottsdale, AZ 85260
Health Canada CTO Registration Certificate Number: 100256

Source Establishment, Processing, Determines Tissue Meets Release Criteria & Distribution: Allotech, LLC
4620 71st Street, Bldg 78-79 Lubbock, TX 79424
Health Canada CTO Registration Certificate number: 100211

Distribution:
Nobel Biocare USA, LLC
1889 West Mission Boulevard
Pomona, CA 91766
Tel.: 1-800-322-5001
info.usa@nobelbiocare.com

Distribution:
Osteogenics Biomedical, Inc.
4620 71st Street, Bldg 78-79
Lubbock, TX 79424
Tel.: 1-806-796-1923
osteogenics.com

CONTACT INFORMATION: Please contact Allotech at 1-888.796.1923 or 1-806.796.1923 to promptly report any unanticipated or adverse events, or should you require further information. Retain this information for your records.

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