



# ZMATRIX™

*porcine peritoneum collagen membrane*

## PERFECTLY **SOFT**


▶ *Zmatrix™ is a natural, native collagen membrane; cross-linking chemicals and agents are unnecessary. Proprietary processing technology allows preservation of collagen as well as extracellular components including laminin, fibronectin, elastin, and glycosaminoglycans.\**

\*Hoganson DM, Owens GE, O'Doherty EM, Bowley CM, Goldman SM, Harilal DO, Neville CM, Kronengold RT, Vacanti JP. Preserved extracellular matrix components and retained biological activity in decellularized porcine mesothelium. *Biomaterials*. 2010, 27: 6934-6940.





- ◀ *Designed to drape without adhering to itself*
- ◀ *Natural peritoneum collagen structure provides for elasticity*
- ◀ *Processed to preserve extracellular components including laminin, fibronectin, elastin, and glycosaminoglycans\**



### Zmatrix™ Porcine Peritoneum Collagen Membrane

15 mm x 20 mm | ZM1520 (1 per box)

20 mm x 30 mm | ZM2030 (1 per box)

30 mm x 40 mm | ZM3040 (1 per box)

- Substantially resorbed within 26 weeks



shown actual size