

10004 - Aatrex 4I

20000 - Aatrex 4L

30000 - DKC60-16 (YGCB)

10007 - Aatrex 4L

3000

55 G

Fille

Inspe

100-817

Bicep II Magnum

Filled '

Not Filled

Cancel

Container ID:

Capacity:

Fill Status:

Product:


EPA#:

Last Activity:

Current Location:



Container Management On The Go



From receiving, filling, or unloading to inspections, cleanouts, or routine maintenance, all activities are recorded in real-time so container information is always up to date.

## **CONTAINER ACTIVITIES**

- Access the app from both iOS and Android phones and tablets
- Customize the container list through settings/filters
- Set the location for containers
- Search or use bar code scanning to find a container
- Perform container activities recording activity dates, times, notes, products, quantities or indicating passed tests, as applicable

Inspect
Fill
Receive from Customer
Cleanout
Information

Maintenance Unload Scheduled Test Remove from Service



Watch the video:

